

PROVAC | Portable Phlegm Suction Machine

Portable absorb phlegm apparatus designed based on developing orientation of similar products at home and abroad is a new generation of oil free lubrication suction device, which is suitable for use by the patient who has difficulty in phlegm removal due to illness, coma and operation, as well as for aspirating such liquid as pus and blood during the clinical practice. It is the commonly applied medical device for use in the emergency room, operation room, and for nursing in sickroom and home health care.

- ❶ With a low flow, medium vacuum and fashionable stylish
- ❷ Adopting oil-free piston pump
- ❸ Adjustable pressure button to select any suitable pressure level
- ❹ Overflow protection can safely and efficiently prevent the liquid from entering into the pump

Main Technical Specifications

1. Power: AC 230V, 50Hz	2. Input current: 0.7A
3. Max negative pressure value : 87 Kpa (0.87 bar)	4. pressure adjustment range: 0.01 MPa ~ Max Pressure Value Limit
5. Pumping rate: ≤ 20 L/MIN	6. Fuse: AC250V, F2.0 AL, Φ5 X 20
7. Reservoir bottle: 1000mL, One Piece	8. Noise: ≤60dB(A)
9. Vacuum gauge: 0~100KPa, 2.5 grade	10. Equipment Type: Medium vacuum / low flow
11. Working mode: Short-time working	12. Electrical requirements: is a class II ,type BF

Note: It is not suitable for use in environment with flammable and explosive gases; Electromagnetic interference source: the minimum separation distance of ≥1.1m when it is close to power lines and electrical equipment with the voltage less than 500V

Accessories list

Suction tubing(length1.3m)	1 PCS	Intermediate tubing (length 0.2m)	2 PCS
Manual	1 PCS	Air filter	3 PCS
Reservoir bottle	1PCS	suction catheter(12F/8F)	2PCS
Catheter connector	1PCS		

